

Superior Thumbkeeper with D-Ring

Maintains the Thumb in Opposition

Application Instructions

1. Loosen straps and slide the thumb into the splint. Secure the thumb with hook and loop closure so there is a snug fit. (Figure 1)




Figure 1

2. The splint should contour the thenar eminence and the hypothenar eminence on the palm of the hand. If it goes beyond this point, trim excess material with a scissors. Secure palmar edge by pulling the strap through the D-ring and securing hook and loop closure. Cut off excess strap. (Figure 2)



Figure 2

3. The thumb splint should not restrict flexion of the IP joint. Trim with scissors if necessary. (Figure 3)



Figure 3

4. There should be a snug fit, but not too tight. Check for red marks once support has been worn. Readjust if necessary.

Note: Right and left splints are not interchangeable.

Washing Instructions: Hand wash with mild soap. Splint may be bent temporarily for cleaning. Towel or air dry.

Note To Physicians and Therapists: The splint may be trimmed with an ordinary pair of scissors. If additional adjustment is needed, spot heat the area with a heat gun, slowly moving it back and forth over the desired area. This will allow you to reshape the splint by applying pressure with your finger.

Reorder No. Left	Reorder No. Right	Size
0814 4952	0814 4962	Small (2 1/4" - 2 3/4")
0814 4953	0814 4963	Medium (2 3/4" - 3 1/4")
0814 4954	0814 4964	Large (3 1/4" - 3 3/4")

Measure width of band at MP joints from radial aspect of index finger to the ulnar aspect of little finger

BIRD & CRONIN INC.

Designers and Manufacturers of Quality Medical Products

1200 Trapp Road • Eagan, Minnesota USA 55121

651/683-1111 • Fax 651/688-9855 • 800/328-1095 • Fax 800/279-7934

e-mail sales@birdcronin.com • www.birdcronin.com

©2000 Bird & Cronin, Inc. All Rights Reserved - Form 4952