

## Maternity-Mate™ Support



Figure 1



Figure 2



Figure 3



Figure 4



Figure 5

**NOT MADE WITH  
natural rubber latex**

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

### Intended for single patient use.



**WARNING:** If you experience any pain or unusual reactions while using this product, consult your healthcare provider immediately.

### APPLICATION INSTRUCTIONS

1. Center padded lift support under the abdomen. Wrap attached strap around waist and secure to support (Figure 1).
2. Adjust pad as low as possible under the abdomen to provide the best support and lift. Adjust pad as needed to avoid interference with legs when walking or sitting (Figure 2).
3. Wrap back support around waist. Grip both ends, pull and secure to padded lift support (Figure 3).
4. To adjust tension and support, release ends of back support, stretch, and reattach (Figure 5).
5. Attach one end of optional upper strap to back support. Softly stretch strap over top of abdomen and secure end on opposite side (Figure 4).

**NOTE:** Do not position the lift support across the front of the abdomen.

**NOTE:** Optional upper strap helps keep support in place and prevent migration.

**NOTE:** Please refer to sew-in label on device for care instructions.

Please visit our website at [www.birdcronin.com](http://www.birdcronin.com) for additional product information or contact a customer care representative at 800-328-1095.

**NOTICE:** Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.