

# Shoulder Abduction Pillow with Harness




Figure 1

## COMPONENTS

A. Shoulder Abduction Pillow  
B. Circumferential Strap  
C. Shoulder Harness Strap  
D. Arm Support Pillow

E. Humeral Cuff  
F. Wrist Cuff  
G. Comfort Foam Pads



Figure 2



Figure 3



Figure 4



Figure 5



Figure 6

**NOT MADE WITH  
natural rubber latex**

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

**This device is intended for single patient use.**



**WARNING:** If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

## APPLICATION INSTRUCTIONS

1. Open package and separate components. (Figure 1)
2. Position the shoulder harness strap (C) on the unaffected shoulder as shown. (Figure 2)
3. Place shoulder abduction pillow (A) under the affected arm. Secure the shoulder harness strap (C) to the front and back of the pillow with the hook tabs provided. (Figure 3)

**NOTE:** Comfort pads may be attached to the front and back of the shoulder harness strap.

4. Attach arm support pillow (D) to the flat, outside edge of the abduction pillow. (Figure 4) Affected arm should rest comfortably on the support.
5. Wrap humeral cuff (E) around bicep of the affected arm and secure to the pillow. (Figure 5)
6. Wrap wrist cuff (F) around wrist of the affected arm and secure to the pillow. (Figure 5)
7. The circumferential strap can be used to secure the pillow in place. Attach the strap to the pillow by wrapping around the patient's mid-section and back onto the pillow. (Figure 6)

**NOTE:** Please refer to sew-in label on device for care instructions.

Please visit our website at [www.birdcronin.com](http://www.birdcronin.com) for additional product information or contact a customer care representative at **800-328-1095**.

**NOTE:** Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.