

NAP Splint - Nighttime Anterior Positioner



Figure 1



Figure 2



Figure 3



Figure 4

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.

WARNING: If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.



APPLICATION INSTRUCTIONS

WARNING: Walking with the NAP splint is NOT recommended. Non-skid tread is present for minimum ambulation.



1. Unfasten the three straps on the splint. (Figure 1)
2. Slide the foot through the elastic piece so the support rests on top of the foot.
3. Pull toes up, wrap and secure the ankle and toe straps tightly. (Figure 2)
4. Wrap the top strap around calf and secure. (Figure 3) Strap should be tight, yet comfortable.
5. Adjust for comfort, not for tightness or excessive stretch.

Please visit our website at www.birdcronin.com for additional product information or contact a customer care representative at **800-328-1095**.

NOTICE: Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.

**NOT MADE WITH
natural rubber latex**