

## B-Cool® Chest & Ab Binder Wrap

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

**Intended for single patient use.**

### Instructions for Use

1. Place plastic bag containing the B-Cool® Chest & Ab Binder Wrap with pre-inserted gel packs into the freezer for at least 4 hours before using.
2. Remove B-Cool® Chest & Ab Binder Wrap from freezer and place over affected area according to physician instructions.
3. Secure the hook and loop closure with enough stretch to provide comfortable support to the affected area. Re-adjust as needed for comfort.
4. Remove B-Cool® Chest & Ab Binder Wrap if it becomes uncomfortable.
5. If pain increases or persists, discontinue use and contact a healthcare provider.

#### CAUTION:

- Never apply a B-Cool® 2.0 gel pack directly onto the skin.
- This product should only be used under the supervision of a medical professional. If pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.
- Supervision is strongly advised when applying this product on children and geriatric patients. Patients with circulatory problems should consult a physician or healthcare professional prior to any treatment.

**NOTE:** Please refer to sew-in label for care instructions.

Please visit our website at [www.birdcronin.com](http://www.birdcronin.com) for additional product and warranty information or contact a customer care representative at 800-328-1095.

**NOTICE:** This device is only one element in the overall treatment program administered by a medical professional. While every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.



## B-Cool® Chest & Ab Binder Wrap

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

**Intended for single patient use.**

### Instructions for Use

1. Place plastic bag containing the B-Cool® Chest & Ab Binder Wrap with pre-inserted gel packs into the freezer for at least 4 hours before using.
2. Remove B-Cool® Chest & Ab Binder Wrap from freezer and place over affected area according to physician instructions.
3. Secure the hook and loop closure with enough stretch to provide comfortable support to the affected area. Re-adjust as needed for comfort.
4. Remove B-Cool® Chest & Ab Binder Wrap if it becomes uncomfortable.
5. If pain increases or persists, discontinue use and contact a healthcare provider.

#### CAUTION:

- Never apply a B-Cool® 2.0 gel pack directly onto the skin.
- This product should only be used under the supervision of a medical professional. If pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.
- Supervision is strongly advised when applying this product on children and geriatric patients. Patients with circulatory problems should consult a physician or healthcare professional prior to any treatment.

**NOTE:** Please refer to sew-in label for care instructions.

Please visit our website at [www.birdcronin.com](http://www.birdcronin.com) for additional product and warranty information or contact a customer care representative at 800-328-1095.

**NOTICE:** This device is only one element in the overall treatment program administered by a medical professional. While every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.



# ***BIRD & CRONIN***<sup>®</sup>

Improving lives with quality medical products

## **B-Cool<sup>®</sup> Chest & Ab Binder Wrap**



PATIENT NAME: \_\_\_\_\_



# ***BIRD & CRONIN***<sup>®</sup>

Improving lives with quality medical products

## **B-Cool<sup>®</sup> Chest & Ab Binder Wrap**



PATIENT NAME: \_\_\_\_\_

